

Determination of Amino Acids in Wine Samples by Precolumn Derivatisation using the FMOc Procedure

Abstract

The application, 'Amino Acid Analysis by Precolumn Derivatisation using a New FMOc Procedure' was used for the amino acid analysis of wine samples, supplied by a local vineyard. This method does not require any solvent extraction or incubation and gives stable single adducts.

Keywords:

Amino Acid, AMINOMATE, Automated Precolumn Derivatisation, FMOc, wine

'...this method does not require any solvent extraction or incubation and gives stable single adducts...'

Sample Preparation

Diluted wine sample (1:10) with derivatisation (diluent) buffer, and filtered.

Conditions

Column: Hypersil Column for Amino Acid Analysis, 150 x 4.6 mm ID

Fluorescence Detector:


Ex: 270 nm
Em: 316 nm

Equilibration: 5 minutes

Flow Rate: 1.00 ml/min

Temperature: 38°C


Injection Vol.: 5 µl



GBC HPLC Instrumentation

LC1150 Quaternary Gradient Pump
LC1120/1150 Online Degasser
LC1250 Fluorescence Detector
LC1650 ACSIS Advanced Autosampler
WinChrom 1-2 Data Management System*

(*plus compatible 486 PC & accessories)



B13
01-0346-00

GBC Scientific Equipment Pty Ltd
A.C.N. 005 472 686
12 Monterey Road, Dandenong, Victoria, 3175, Australia
Phone: (03) 9213 3666 Fax: (03) 9213 3677

All rights reserved
GBC publication number 01-0346-00
September, 1995