

Determination of Ammonia in Water by Pre-Column Derivatisation

Abstract

A semi-automated method is described for the analysis of ammonia in water by pre-column derivatisation with FMOCl (Fluorenylmethyl chloroformate), followed by reversed phase chromatography with UV detection. This method is highly selective and sensitivity is at the 0.1 ppm level. The derivatisation procedure is straightforward and no costly reagent is required.

Methods for the analysis of ammonia in water are numerous, ranging from relatively simple colorimetric methods¹ to sophisticated microcomputer-controlled flow-analysis system². Among various procedures, HPLC offers the versatility of modifying a standard LC procedure to suit the requirements of particular analyses. These requirements include criteria such as sensitivity, selectivity and degree of instrument automation.

FMOCl has been employed for the derivatisation of amino acids³. This derivatisation procedure does not require any solvent extraction and gives stable derivatives. An automated analysis system based on this chemistry has also been reported⁴. By employing similar chemistry, a method has been developed for the selective analysis of ammonia in water by pre-column derivatisation, followed by reverse-phase HPLC with binary gradient and UV detection. In addition, the procedure offers the possibility of further automation, including sample preparation by customisation of the GBC WinChrom Data Station and LC1610 Advanced Autosampler.

Keywords:

Ammonia, Water, Pre-column Derivatisation, FMOCl




Figure 1 Ammonia Standard (5 ppm)

Conditions

Column: Spherisorb S5 ODS₂,
150 x 4.6 mm ID

Mobile Phase: Solvent A - Water, Solvent B -
methanol (Helium Sparging)

Gradient Program:

Time (min)	A (%)	B (%)
0.0	60.0	40.0
1.0	60.0	40.0
25.0	35.0	65.0
26.0	0.0	100.0
30.0	0.0	100.0
31.0	60.0	40.0

Injection Vol: 20 µl
Flow Rate: 1.50 ml/min
Temperature: Ambient
Detection: UV at 263 nm

'...FMOCl
has been
employed for
the
derivatisation of
amino acids...'



*'...by employing
similar
chemistry, a
method has
been developed
for the selective
analysis of
ammonia...'*

Derivatisation

Derivatisation Buffer:

0.25 M Sodium Borate, pH 9.6

Derivatisation Reagent:

FMOC-Cl (25 mg/ml) in Acetonitrile (prepare daily)

Derivatisation Procedure:

A sample of water is collected under clean conditions and filtered through a 0.2 μ filter. A 100 μ l aliquot is added to a 1.5 ml Eppendorf tube. To this is added 450 μ l of borate buffer and 450 μ l of acetonitrile, followed by 50 μ l of the derivatisation reagent.

The vial is capped and allowed to react at room temperature for 5 minutes. At the end of the reaction period, 10 μ l of glacial acetic acid is added to the mixture to quench the reaction. The mixture is recapped and mixed thoroughly. A 20 μ l aliquot of the mixture is subjected to HPLC analysis.

GBC HPLC Instrumentation

LC1110 Dual Piston HPLC Pump (x 2)

LC1200 Variable Wavelength UV/Vis

Detector

LC1431 System Organiser

LC1650 Advanced Autosampler

WinChrom Chromatography Data

Management System





Figure 2 Analysis of Ammonia in Hail Water

References

1. M.A. Brzezinski, *Mar. Chem.*, 20(3), 1987, 277.
2. H. Hara, A. Motoike and S. Okazaki, *Anal. Chem.*, 59(15), 1987, 1995.
3. P.A. Haynes, D. Sheumack, J. Kibby and J.W. Redmond, *J. Chromatogr.*, 540, (1991), 177.
4. P.A. Haynes, D. Sheumack, L.G. Greig, J. Kibby and J.W. Redmond, in press.



E3
01-0349-00

GBC Scientific Equipment Pty Ltd

A.C.N. 005 472 686

12 Monterey Road, Dandenong, Victoria, 3175, Australia

Phone: (03) 9213 3666 Fax: (03) 9213 3677

All rights reserved

GBC publication number 01-0349-00

September, 1995