


Benzidine Separation Using Glassy Carbon Electrode (3 mm) with ppb Sensitivity

'...the sample preparation in this case consists of filtering of the sample, followed by injection...'

Peak 1: Benzidine (10 ppb)
Peak 2: 3,3'-Dichlorobenzidine (25 ppb)



Conditions

Column: Spherisorb 5 μ OD S2 (C18), 150 x 4.6 mm ID
Temperature: Ambient
Mobile Phase: 50% 0.1 M Sodium Acetate (pH 4.7) with 50% Acetonitrile
Flow Rate: 1.5 ml/min
Electrode: 3 mm Glassy Carbon
Potential: 800 mV (Ag/Ag/Cl)

Key Features

- Sensitivity
- Cost savings
- Selectivity
- New patented technology
- Limited technical expertise

Relevant Industries

- Commercial environmental laboratories
- EPA laboratories
- CLP
- Industrial monitoring laboratories
- Water quality laboratories
- Water treatment laboratories
- Landfill monitoring laboratories


This application illustrates the use of the patented impinging wall jet electrochemical detector. Most environmental applications require extractions and derivative formation to ensure sensitivity and selectivity. This detector removes labour intensive sample preparation required for other electrochemical detectors. This can mean up to a days labour saving. The sample preparation in this case consists of filtering of the sample, followed by injection. When sample throughput is required, this application can save the commercial clinical laboratory hundreds of dollars per run.



E11
01-0357-00

GBC HPLC Instrumentation

LC1110 Dual Piston HPLC Pump
LC1260 Electrochemical Detector
LC1445 System Organiser
LC1650 Advanced Autosampler
WinChrom Chromatography Data
Management System



E11
01-0357-00

GBC Scientific Equipment Pty Ltd
A.C.N. 005 472 686
12 Monterey Road, Dandenong, Victoria, 3175, Australia
Phone: (03) 9213 3666 Fax: (03) 9213 3677

All rights reserved
GBC publication number 01-0357-00
September, 1995