



GBC Scientific Equipment

Pre-Owned FS3000 Autosampler

\$2,500.00



FS3000 Autosampler

Automatic sample changer for flame analysis

- Carousel holds 10 standards and 60 samples.
- Separate positions for rinse, blank and rescale standard.
- Programmable measurement time, delay time, rinse time, rescale rate, recalibration rate, number of standards, number of replicates, first and last sample numbers.

For use with GBC Flame AAS instruments ONLY



GBC Scientific Equipment (USA) LLC
151A N. State Street, P.O. Box 339
Hampshire, IL 60140
Toll-Free: 800-445-1902 Fax: 847-683-9871
Email: sales@gbcscientific.com
www.gbcscientific.com